

Analyze D₂O Losses from the D₂O loaded AgZ provided in FY15 By ORNL

Fuel Cycle Research & Development

Milestone M3FT-16SN030107041

DOE/Ne-Fuel Cycle R&D
Materials Recovery and Waste Form Development

Prepared for
US Department of Energy
Office of Nuclear Energy –
Materials Recovery and Waste Form Development
Todd Alam, Tina M. Nenoff (Corresponding
Author), Russel Jarek, Curtis D. Mowry
Sandia National Laboratories
April 20, 2016
FCRD-MRWFD-2016-000259



DISCLAIMER

This information was prepared as an account of work sponsored by an agency of the U.S. Government. Neither the U.S. Government nor any agency thereof, nor any of their employees, makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness, of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. References herein to any specific commercial product, process, or service by trade name, trade mark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the U.S. Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the U.S. Government or any agency thereof.

Develop Plan for Analysis of the Effluent from GCM Production

Fuel Cycle Research & Development

Prepared for
U.S. Department of Energy
Materials Recovery and Waste Form Development

Todd Alam, Tina M. Nenoff, Russel Jarek, Curtis Mowry

Sandia National Laboratories

April 20, 2016

FCRD-MRWFD-2016-000259

SAND 2016-XXXX

Plan:

This milestone is focused on utilizing NMR and various gas chromatograph set ups (thermal conductivity and/ionizing detection) at SNL, to analyze the relative ratio of water and deuterated water in Ag-MOR samples of natural zeolite sent to SNL in FY15 by ORNL.

Samples Studied:

ORNL D₂O AgI-MOR sample(s) were I₂ and D₂O loaded, then additional treatments of:

Ag[°]Z-I-D₂O SHB #1: 24 hr dry air purge

Ag[°]Z-I-D₂O SHB #2: 4 hr dry air purge, moist air stream loading, 24 hr dry air purge

Analysis Methods:

1) The solid state MAS NMR Bruker Avance – III 600 NMR instrument using a 2.5 mm rotor spinning at 30 kHz for proton, and a 4 mm rotor spinning at 12 kHz for deuterium characterization. The ¹H NMR chemical shifts were referenced to the secondary external standard adamantane $\delta = 1.63$ ppm with respect to TMS = 0.0 ppm. Deconvolutions were obtained using the DMFIT software.

2) Isotope measurements via Mass Spec (Finnegan MAT 271 magnetic sector mass spec). Instrument has a limitation of inherent surface composition of –OH's; even at ~100 milliTorr, the -OH levels were swamping the very low level of D₂O (if any) in the sample.

3) Isotope measurements via pulsed discharge ionization detector (PDID; Valco Instruments Co., VICI AG International), to separate D₂O and H₂O by Gas Chromatography.

Results:

NMR: In the two materials, the ¹H (proton) NMR signal is dominated by the broad water resonance near ~4.5 ppm. This resonance reflects waters that are not strongly hydrogen bonded to the surface or other functional groups, but are surface adsorbed resulting in short T2 values (see below) and broader lines. It is actually heterogeneous in nature and can be considered the overlap of several different water environments and slightly different bonding environments within the pores (See integrations in Figures 5 – 9). These results also show that the SNL D2O-X samples have the lowest concentration of water and other proton species in comparison to the other samples. The ORNL AGI-MOR sample has a the sharp resonance at ~ 9 ppm, which reflects a proton (OH of H₂O) in a relatively strong hydrogen bonding environment.

For the deuterium analysis (²H) of the two ORNL samples: The chemical shift is consistent with water in these materials. The presence of multiple sidebands demonstrates that the water is rigid and not isotopically moving. The goal we were tasked with was to determine the relative concentration of ²H remaining in these materials. The presence of deuterium is small in either of these samples, only a ratio between them is measurable. Since the *total integration is related to the concentration, the ratio of the D₂O-2/D₂O-1 in these samples was found to be 13% or about a 7.7 fold reduction in the ²H concentration following the additional air flow treatment.*

Isotope Measurement and PDID: Neither of these analyses were able to produce results to determine concentrations of D₂O in the samples. In fact, the only conclusive result here is that the quantity of deuterium was too small to detect.

ACKNOWLEDGEMENTS

The authors thank S. Bruffey and R. Jubin (ORNL) for supplying the D₂O-AgI-Mordenite samples.

We acknowledge funding from the DOE/NE Fuel Cycle R&D – Materials Recovery and Waste Form Development (MRWFD).

Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.

DISTRIBUTION LIST

- 1 MS1415 Tina M. Nenoff, SNL, 1100
- 1 MS0886 Curt D. Mowry, SNL, 1852
- 1 MS0889 Katheryn Helean, SNL, 1852
- 1 MS1427 Grant Heffelfinger, SNL, 1100
- 1 MS0899 Technical Library, SNL, 9536 (electronic copy)

- 1 James Bresee
Office of the Asst. Sec. for Nuclear Energy
Global Nuclear Energy Partnership (NE-2.4)
U.S. Department of Energy
1000 Independence Ave. SW
Washington, DC, 20585

- 1 Kimberly Gray
Office of the Asst. Sec. for Nuclear Energy
Global Nuclear Energy Partnership (NE-2.4)
U.S. Department of Energy
1000 Independence Ave. SW
Washington, DC, 20585

- 1 Terry A Todd
Idaho National Laboratory
PO Box 1625
Idaho Falls, ID 83415

- 1 John Vienna
Pacific Northwest National Laboratory
902 Battelle Blvd.
PO Box 999
Richland, WA 99352

- 1 Robert Jubin
Oak Ridge National Laboratory
1 Bethel Valley Road
P.O. Box 2008, MS 6243
Oak Ridge, TN 37831-6243

- 1 TIO Document Management
GNEPTIODMS@inl.gov