
## PEEM-iv vs hu-scan: case study for SAND2017-2314R graphene folds on SiO2

- Determining the work function and the highest occupied state at Γ-point from PEEM-iv
  - Highest occupied state at Γ-point corresponds to the Fermi level for graphene due to its metallic character
  - PEEM-iv measured using  $\lambda = 190$  nm,  $4^{th}$ analyzer slit (~200 meV resolution)

Determining ionization energy from the photoemission threshold from photon energy scan (hu scan)

 $\phi$  ~4.60 eV: sample is

Ionization energy corresponds to the work function for metallic material




## PEEM-iv vs hu-scan: case study for 2ML graphene island on SiO2

- The vacuum level cut-off data shows ~80meV and ~70meV difference between 1GL vs 2GL and 2GL vs multi-layer
  - These variations are in a good agreement with Fig. 4 in the PRB
- The upper part of the 2GL island is the slight variation (~30meV), which I interpret as a twist angle dependence of the work function in 2GL

0.05 -0.00 --0.05-0.10 -0.15 -0.1 0.0 (eV) -0.1-0.230µm -0.05 -0.10~70meV -0.15

The work function of 1GL, ~4.43eV is in good agreement with the calculated work function of 1GL with mid 10^12 cm-2 doping shown in Fig. 5 in the PRB.



## • User facility:

 Center for Integrated Nanotechnologies (CINT), SNL, supported by, an Office of Science User Facility operated for the U.S. DOE Office of Science by Sandia National Laboratories (Contract DE-AC04-94AL85000) and Los Alamos National Laboratory (Contract DE-AC52-06NA25396).



Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.